



Early View

Original research article

Rapid, simplified whole blood-based multiparameter assay to quantify and phenotype SARS-CoV-2 specific T cells

Catherine Riou, Georgia Schäfer, Elsa du Bruyn, Rene T. Goliath, Cari Stek, Huihui Mou, Deli Hung, Katalin A. Wilkinson, Robert J. Wilkinson

Please cite this article as: Riou C, Schäfer G, du Bruyn E, *et al.* Rapid, simplified whole blood-based multiparameter assay to quantify and phenotype SARS-CoV-2 specific T cells. *Eur Respir J* 2021; in press (<https://doi.org/10.1183/13993003.00285-2021>).

This manuscript has recently been accepted for publication in the *European Respiratory Journal*. It is published here in its accepted form prior to copyediting and typesetting by our production team. After these production processes are complete and the authors have approved the resulting proofs, the article will move to the latest issue of the ERJ online.

Copyright ©The authors 2021. This version is distributed under the terms of the Creative Commons Attribution Non-Commercial Licence 4.0. For commercial reproduction rights and permissions contact permissions@ersnet.org

Rapid, simplified whole blood-based multiparameter assay to quantify and phenotype SARS-CoV-2 specific T cells

Catherine Riou^{1,2,*}, Georgia Schäfer^{1,3,4}, Elsa du Bruyn^{1,5}, Rene T. Goliath¹, Cari Stek^{1,5,6}, Huihui Mou⁷, Deli Hung⁸, Katalin A. Wilkinson^{1,5,9}, Robert J. Wilkinson^{1,5,6,9}

¹ Wellcome Centre for Infectious Disease Research in Africa and Institute of Infectious Disease and Molecular Medicine, University of Cape Town, Observatory 7925, South Africa

² Division of Immunology, Department of Pathology, University of Cape Town, Observatory 7925, South Africa

³ International Centre for Genetic Engineering and Biotechnology (ICGEB) Cape Town, Observatory 7925, South Africa

⁴ Division of Medical Biochemistry and Structural Biology, Department of Integrative Biomedical Sciences, University of Cape Town, Observatory 7925, South Africa

⁵ Department of Medicine, University of Cape Town, Observatory 7925, South Africa

⁶ Department of Infectious Diseases, Imperial College London, W12 0NN, UK

⁷ Department of Immunology and Microbiology, The Scripps Research Institute, Jupiter, FL 33458, USA

⁸ Department of Immunology and Microbiology, The Scripps Research Institute, La Jolla, CA 92037, USA

⁹ The Francis Crick Institute, London, NW1 1AT, UK

* **Corresponding Author:** Catherine Riou, University of Cape Town, IDM, Anzio Road, Observatory, 7925, Cape Town, South Africa. E-mail: cr.riou@uct.ac.za

Key words: COVID-19, T cells, Whole blood Assay

summary: In this proof-of-concept study, we show that SARS-CoV-2 T cell responses are easily detectable using a rapid whole blood assay requiring minimal blood volume. Such assay represents a suitable tool to monitor adaptive immunity in vaccine trials.

ABSTRACT

Rapid tests to evaluate SARS-CoV-2-specific T cell responses are urgently needed to decipher protective immunity and aid monitoring vaccine-induced immunity. Using a rapid whole blood assay requiring minimal amount of blood, we measured qualitatively and quantitatively SARS-CoV-2-specific CD4 T cell responses in 31 healthcare workers, using flow cytometry. 100% of COVID-19 convalescent participants displayed a detectable SARS-CoV-2-specific CD4 T cell response. SARS-CoV-2-responding cells were also detected in 40.9% of participants with no COVID-19-associated symptoms or who tested PCR negative. Phenotypic assessment indicated that, in COVID-19 convalescent participants, SARS-CoV-2 CD4 responses displayed an early differentiated memory phenotype with limited capacity to produce IFN γ . Conversely, in participants with no reported symptoms, SARS-CoV-2 CD4 responses were enriched in late differentiated cells, co-expressing IFN γ and TNF α and also Granzyme B. This proof-of-concept study presents a scalable alternative to PBMC-based assays to enumerate and phenotype SARS-CoV-2-responding T cells, thus representing a practical tool to monitor adaptive immunity due to natural infection or vaccine trials.

INTRODUCTION

The outbreak of SARS-CoV-2 infection (causing the disease known as COVID-19) that first emerged in Wuhan, China in December 2019, was declared a global pandemic on 12th March 2020, and is affecting all countries of the world, including those of Africa [1]. There is an urgent need to understand better the clinical manifestations and the pathogenesis of SARS-CoV-2 in order to develop relevant tools, including diagnostic tests, treatments and vaccines, to stop the spread of disease as well as strategies to best manage this disease in all population groups.

The clinical spectrum of COVID-19 is very wide, from asymptomatic through mild flu-like symptoms to severe pneumonia and death. Understanding what constitutes immune protection against SARS-CoV-2 is key to predicting long-term immunity and to inform vaccine design. While much emphasis has been placed on the B cell and antibody response, it is not yet clear what type of immune response confers protection to SARS-CoV-2 [2]. Several studies suggest that the T cell response may play an important role in SARS-CoV-2 pathogenesis and reports indicating that patients lacking B cells can recover from SARS-CoV-2 infection further highlight the likely importance of T cell immunity [3, 4]. Additionally, accumulating evidence indicates that the presence of pre-existing, cross-reactive memory T cells specific for common cold coronaviruses may affect susceptibility to SARS-CoV-2 infection and partially explain the markedly divergent clinical manifestations of COVID-19 [5-8].

Efficient, sensitive and simple assessment of human T-cell immunity remains a challenge. Most commonly used T cell assays necessitate the isolation of peripheral blood mononuclear cells (PBMC), requiring significant amount of blood. Therefore, whole blood assays could be more advantageous than PBMC-based methods, by significantly reducing blood volume (~ 1 ml), making them more applicable to paediatric populations. Moreover, such assays are rapid, as they

don't require cell separation, and preserve the physiologic cellular and soluble environments, mimicking better human blood condition. Here we report a rapid (~ 7 hrs) whole blood-based detection method of SARS-CoV-2-specific T cell responses with simple steps that could be adapted to settings of limited resources. This rapidly applicable assay could represent an easily standardizable tool to assess SARS-CoV-2-specific adaptive immunity to monitor T cell responses in vaccine trials, gain insight into what constitutes a protective response, or define the prevalence of SARS-CoV-2 T cell responders population wide.

METHODS

Study population

The population studied consisted of healthcare workers (HCW, n=31, 29% male) recruited between July and September 2020, from Groote Schuur Hospital in Cape Town, the hardest hit region of the initial COVID-19 epidemic in South Africa [9]. Participants were classified according to i) reporting of COVID-19-associated symptoms, ii) whether a SARS-CoV-2 RNA PCR from a nasal or pharyngeal swab was performed and iii) SARS-CoV-2 PCR results. Based on these criteria, participants were subdivided into three groups: 1) persons with no COVID-19-associated symptoms (n=15); 2) persons who reported symptoms but tested SARS-CoV-2 PCR negative (n = 7); and 3) persons who had COVID-19-associated symptoms and tested SARS-CoV-2 PCR positive (n = 9) (**Table 1**). In all groups, the exposure to COVID-19 patients was comparable (86.8 to 100%). All participants with PCR confirmed COVID-19 had mild symptoms and did not require hospitalization. Blood samples were obtained a median of 7.3 weeks post-SARS-CoV-2 testing in persons with a negative test results and 4.7 weeks in those with a confirmed positive result (p=0.09). All participants were symptom-free at the time of sampling. In four SARS-CoV-2 PCR positive participants, a second sample was obtained 1 month later. The University of Cape Town's Faculty of Health Sciences Human Research Ethics Committee

approved the study (HREC: 207/2020) and written informed consent was obtained from all participants.

Measurement of SARS-CoV-2 nucleocapsid-specific IgG in plasma

The measurement of SARS-CoV-2 specific antibodies was performed using the Roche Elecsys[®] Anti-SARS-CoV-2 immunoassay (Roche Diagnostics, Basel, Switzerland). This semi-quantitative electro-chemiluminescent immunoassay measures SARS-CoV-2 nucleocapsid-specific IgG. The assay was performed by the South African National Health Laboratory Service (NHLS) and interpreted according to manufacturers' instructions (Roche: V 1.0 2020-05). Results are reported as numeric values in form of a cut-off index (COI; signal sample/cut-off), where a COI < 1.0 corresponds to non-reactive plasma and COI ≥ 1.0 to reactive plasma. At 14 days post-SARS-CoV-2 PCR confirmation, the sensitivity and specificity of the Elecsys[®] Anti-SARS-CoV-2 immunoassay has been reported as 99.5% (95% CI, 97.0 to 100.0%) and 99.80% (95% CI, 99.69 to 99.88%), respectively [10-12].

SARS-CoV-2 pseudovirus neutralisation assay

Patient plasma was evaluated for SARS-CoV-2 neutralisation activity using a SARS-CoV-2 pseudovirus infection assay. Single-cycle infectious SARS-CoV-2 pseudovirions based on the HIV backbone expressing the SARS-CoV-2 Spike (S) protein and a firefly luciferase reporter were produced in HEK-293TT cells [13, 14] by co-transfection of plasmids pNL4-3.Luc.R-E- (NIH AIDS Reagent Program (#3418), Germantown, MD, USA) and pcDNA3.3-SARS-CoV-2-spike Δ18 [15]. HEK-293TT cell culture supernatants containing the virions were harvested 3 days post transfection and incubated with heat-inactivated patient plasma at 5-fold serial dilutions for 60 min at 37°C. Plasma/pseudovirus mixtures were then used for transfection of HEK-293T cells stably expressing the ACE2 receptor [16]. Cells were lysed 3 days post infection using the Promega cell culture lysis reagent (Promega Biosciences Inc., San Luis Obispo, CA, USA) and

assessed for luciferase activity using a GloMax[®] Explorer Multimode Microplate Reader (Promega Biosciences) together with the Luciferase assay system (Promega Biosciences).

Whole blood-based T cell detection assay

Blood was collected in sodium heparin tubes and processed within 3 h of collection. The whole blood-based SARS-CoV-2-specific T cell detection assay was adapted from a previously reported whole blood intracellular cytokine detection assay designed to quantitate *Mycobacterium tuberculosis* specific T cells in small volumes of blood [17]. However, significant modifications have been made, including a reduced incubation time, the usage of a fixation buffer allowing the simultaneous lysis of red blood cells to streamline processing time, leading to faster acquisition of results [18]. Here, we adapted this assay to detect SARS-CoV-2 specific T cells using synthetic SARS-CoV-2 PepTivator peptides (Miltenyi Biotec, Surrey, UK), consisting of 15-mer sequences with 11 amino acid overlap covering the immunodominant parts of the spike (S) protein, and the complete sequence of the nucleocapsid (N) and membrane (M) proteins. All peptides were combined in a single pool and used at a final concentration of 1 µg/ml. The workflow of the assay is presented in **Figure 1**. Briefly, 400 µl whole blood was stimulated with the SARS-CoV-2 S, N and M protein peptide pool at 37°C for 5 hrs in the presence of the co-stimulatory antibodies against CD28 and CD49d (1µg/ml each; BD Biosciences, San Jose, CA, USA) and Brefeldin-A (10µg/ml, Sigma-Aldrich, St Louis, MO, USA). Unstimulated blood was incubated with co-stimulatory antibodies, Brefeldin-A and an equimolar amount of DMSO. Red blood cell lysis and white cell fixation was then performed as a single step using a Transcription Factor Fixation buffer (eBioscience, San Diego, CA, USA) for 20 minutes. At this stage cells were cryopreserved in freezing media (50% foetal bovine serum, 40% RPMI and 10% dimethyl sulfoxide) and stored in liquid nitrogen until batched analysis.

Cell staining was performed on cryopreserved cells that were thawed, washed and permeabilised with a Transcription Factor perm/wash buffer (eBioscience). Cells were then stained at room temperature for 30 min with antibodies for CD3 BV650, CD4 BV785, CD8 BV510, CD45RA Alexa 488, CD27 PE-Cy5, CD38 APC, HLA-DR BV605, Ki67 PerCP-Cy5.5, PD-1 PE, Granzyme B (GrB) BV421, IFN γ BV711, TNF α PE-Cy7 and IL-2 PE/Dazzle 594, as detailed in **Supplementary Table 1**. Samples were acquired on a BD LSR-II and analysed using FlowJo (v9.9.6, FlowJo LCC, Ashland, OR, USA). A positive response was defined as any cytokine response that was at least twice the background of unstimulated cells. To define the phenotype of SARS-CoV-2-specific CD4 T cells, a cut-off of 30 events was used. The gating strategy is provided in **Supplementary Figure 1**.

Statistical analyses

Graphical representations were performed in Prism (v8.4.3; GraphPad Software Inc, San Diego, CA, USA) and JMP (v14.0.0; SAS Institute, Cary, NC, USA). Statistical tests were performed in Prism. Non-parametric tests were used for all comparisons. The Kruskal-Wallis test with Dunn's multiple comparison test was used for multiple comparisons and the Mann-Whitney and Wilcoxon matched pairs test for unmatched and paired samples, respectively.

RESULTS

Serological assessment of SARS-CoV-2 sensitization.

Even though RT-PCR is the most specific technique to detect acute SARS-CoV-2 infection, the positivity rate drops rapidly as soon as 10 days post-symptom onset, particularly in individuals with mild forms of COVID-19 [19, 20]. Hence, serology assays provide an important complement to RNA testing to identify individuals who have been sensitized by SARS-CoV-2. Thus, to assess potential SARS-CoV-2 sensitization in participants who did not have a SARS-

CoV-2 PCR test performed or tested PCR negative, the presence of SARS-CoV-2-specific antibodies (e.g. SARS-CoV-2 nucleocapsid-specific IgG) was measured and a SARS-CoV-2 pseudovirus neutralisation assay was also performed in all participants. While all participants with a positive SARS-CoV-2 PCR test exhibited *in vitro* anti-SARS-CoV-2 neutralizing activity and 8/9 (88.8%) had detectable SARS-CoV-2 nucleocapsid-specific IgG, none of the participants who tested negative or did not undergo PCR testing were positive for SARS-CoV-2 nucleocapsid-specific IgG (**Figure 2A**) or displayed robust *in vitro* anti-SARS-CoV-2 activity (**Figure 2B**). Overall, these results confirmed that PCR positive participants had been infected with SARS-CoV-2 and mounted an immune response to the virus. For the other participants, it is difficult to ascertain their SARS-CoV-2 infection status. On one hand, the absence of antibody response and neutralization activity may suggest that they have not been SARS-CoV-2 infected. But on the other hand, due to their high level of SARS-CoV-2 exposure, we cannot fully exclude SARS-CoV-2 infection, as while uncommon, COVID-19 confirmed patients with negative antibody results have been reported [21].

Magnitude and functional profile of SARS-CoV-2-responding CD4 T cells.

Amongst the 31 participants tested, 58% (n=18) had a detectable SARS-CoV-2-specific CD4 T cell response (producing any of the measured cytokines, IFN γ , TNF α or IL-2) using the described whole blood assay (**Figure 3A**). We then defined the magnitude and phenotype of SARS-CoV-2-specific CD4 responses according to participants' clinical characteristics (**Figure 3B&C**). All HCW with self-reported symptoms and a SARS-CoV-2 PCR positive test (n=9) were found to have SARS-CoV-2-specific CD4 T cells (median frequency: 0.47%, IQR: 0.28-0.65) (**Figure 3C**). Of note, the frequency of SARS-CoV-2 specific CD4 T cells did not associate with the magnitude of SARS-CoV-2 nucleocapsid-specific IgG (p = 0.64, r = -0.18) or plasma neutralizing activity (p = 0.46, r = 0.28) (data not shown). Interestingly, in participants with no detectable SARS-CoV-2 nucleocapsid-specific IgG or neutralizing antibody activity, of those with self-reported symptoms, 3/7 (43%) had detectable SARS-CoV-2-responding CD4 T cells

(median of responders: 0.02%, IQR: 0.019-0.03) and 6/15 (40%) of those with no self-reported symptoms had detectable SARS-CoV-2-responding CD4 T cells (median of responders: 0.15%, IQR: 0.05-1.47). Moreover, a SARS-CoV-2-specific CD8 T cell response was also observed in 6 out of the 9 SARS-CoV-2 PCR positive participants (66.7%). In participants reporting no symptoms (n=15), only one person had a detectable CD8 response, and none of the PCR negative participants (n=7) exhibited a detectable SARS-CoV-2-specific CD8 T cell response (data not shown).

Next, we defined the polyfunctional profile of SARS-CoV-2-responding CD4 T cells, based on their capacity to co-express IL-2, IFN γ or TNF α (**Figure 3D**). The overall functional profile of SARS-CoV-2-specific cells in PCR positive participants was distinct from participants with no reported symptoms (p=0.0002). Indeed, the SARS-CoV-2-specific CD4 response in PCR+ participants was characterized by limited expression of IFN γ and was enriched in cells co-expressing IL-2 and TNF α . On the contrary, in participants reporting no symptoms, most SARS-responding CD4 cells were distributed between triple functional cells (IL-2+IFN γ +TNF α +) and cells co-producing IFN γ and TNF α . Overall, in this assay TNF α was the predominant cytokine produced, with its production being significantly higher compared to IL-2 (p = 0.0026) and IFN γ (p < 0.0001) (**Supplementary Figure 2**).

Phenotypic assessment of SARS-CoV-2-responding CD4 T cells

While the proposed assay can be performed using a limited antibody panel to identify the frequency of SARS-CoV-2-responding T cells, by solely measuring cytokine production, the use of a more extensive antibody panel also permits definition of the phenotypic profile of these cells. To this end, we included additional markers to assess the memory differentiation (CD27, CD45RA), cytotoxic potential (GrB) and activation profile (HLA-DR, CD38, Ki67, PD-1) of SARS-CoV-2-responding CD4 T cells.

Figure 4A&B shows that in all individuals with symptoms, SARS-CoV-2-specific CD4 T cells displayed almost exclusively (median: 97.7%) an early differentiated memory phenotype (ED: CD45RA-CD27+). On the contrary, in 4/6 HCW with no symptoms, SARS-CoV-2-specific CD4 T cells exhibited a predominant late differentiated phenotype (LD: CD45RA-CD27-).

While the role of cytotoxic CD4 T cells is still unclear, the presence of these cells has been described in several viral infections [22]. We thus measured GrB expression in SARS-CoV-2-responding CD4 T cells. While GrB was barely detectable in SARS-CoV-2-specific CD4 T cells in PCR positive participants, elevated GrB expression was observed in 4 out of 6 participants without symptoms (**Figure 4C&D**). Moreover, the proportion of ED SARS-CoV-2-responding CD4 T cells inversely associated with GrB expression ($p = 0.002$, $r = -0.71$) and the proportion of IFN γ and TNF α dual producing cells ($p < 0.0001$, $r = -0.84$) (**Figure 4E**). The phenotypic characteristics of GrB expressing SARS-CoV-2-responding CD4 T cells (i.e. highly differentiated and producing IFN γ and TNF α) are in agreement with a recent report describing CMV-specific CD4 CTL T cells [23].

To determine if the overall phenotypic profile of SARS-CoV-2-responding CD4 T cells allowed discrimination of participants based on their clinical characteristics, we performed a hierarchical clustering analysis (**Figure 5A**) and a principal component analysis (PCA) (**Figure 5B**), including four parameters (e.g. proportion of IFN γ +TNF α +IL2+ and IFN γ -TNF α +IL2+ cells, proportion of ED, and GrB expression). Both analyses show that participants who reported no symptoms and were negative for SARS-CoV-2 Abs separated clearly from PCR positive participants. Conversely, PCR and SARS-CoV-2 antibody negative HCW reporting symptoms could not be separated from PCR positive participants. The loading plot shows that GrB expression and the proportion of triple positive cells were the main drivers permitting the segregation of participants who reported no symptoms and were negative for SARS-CoV-2 Abs (**Figure 5C**)

Lastly, as the expression of activation markers (HLA-DR, CD38, PD-1 or Ki67) on antigen-specific T cells are indicative of active infection [24, 25], we defined the expression of these markers on SARS-CoV-2-responding CD4 T cells. No significant difference was observed in any of the measured markers amongst the different groups, notwithstanding a few outliers observed in the PCR positive participant group (**Figure 6A**). Interestingly, regardless of the participants' clinical characteristics, PD-1 expression was highly expressed on SARS-CoV-2-responding CD4 T cells, We then defined the relationship between the activation profile of SARS-CoV-2-specific CD4 T cells and time post positive PCR test. **Figure 6B** shows that in PCR positive participants, the expression of CD38 associated strongly with the time post PCR testing ($p < 0.0001$, $r = -0.91$), where CD38 expression decreased sharply in the first 2 to 3 weeks post testing, suggesting a rapid clearance of the pathogen. Similar associations were observed for all tested activation markers, with CD38 expression showing the strongest correlation (**Figure 6C**).

DISCUSSION

In this proof-of-concept analysis assessing the use of a simple whole blood assay to measure SARS-CoV-2-specific T cell responses in a small health care worker cohort, we show that SARS-CoV-2-specific CD4 T cells (expressing IFN γ , TNF α or IL-2) were easily detectable in all tested convalescent COVID-19 participants. However, in participants who did not experience any COVID-19-related symptoms or tested SARS-CoV-2 PCR negative (despite reporting symptoms), 9/22 (40.9%) also exhibited a detectable SARS-CoV-2-specific CD4 T cell response. Nevertheless, the median frequency of SARS-CoV-2-specific CD4 T cells in the latter groups was ~ 5-fold lower compared to SARS-CoV-2 PCR positive participants. Importantly, unlike convalescent COVID-19 participants, none of participants without symptoms or who tested PCR negative had detectable SARS-CoV-2 nucleocapsid-specific IgG or displayed robust neutralizing activity. Several hypotheses could explain the detection of SARS-CoV-2 responding cells in the absence of detectable antibody response. As the participants included in this study are most likely highly exposed to SARS-COV-2 due to their occupation, one could argue that 1) SARS-CoV-2 CD4 T cell response was generated during an asymptomatic COVID-19 episode, that did not result in a detectable antibody response; 2) reflect recent exposure to SARS-CoV-2; or 3) correspond to SARS-CoV-2 cross-reactive memory CD4 T cells. A number of studies have demonstrated the presence of SARS-CoV-2-reactive CD4 T cells in 40 to 60% of SARS-CoV-2-unexposed individuals in different populations around the world [26-29]. Additionally, it is important to point out that all these studies were performed using PBMC and it is encouraging to find that a whole blood-based assay using limited amount of blood (< 1ml) yields comparable results.

Further analyses of the polyfunctional and phenotypic profile of SARS-CoV-2-specific CD4 T cells revealed that regardless of the clinical characteristics of the participants, the most prevalent cytokine detected in response to SARS-CoV-2 peptides was TNF α . These observations

are in accordance with other studies pointing out an impairment of SARS-CoV-2 adaptive immune responses, characterized by low level of type I and type II interferon [30, 31]. However, it is also possible that the limited IFN γ production, we observed in this study, is related to the short stimulation time (e.g. 5 hrs) used in the assay. Nevertheless, this suggests that TNF α could be a more reliable target than IFN- γ , to detect SARS-CoV-2-specific CD4 T cell responses. Additionally, SARS-CoV-2-responding CD4 T cells in all participants were characterized by elevated of PD-1 expression. Such profiles have been observed in acute SARS-CoV-2 infection and is likely driven by ongoing viral replication [27, 32, 33]. Expression of PD-1 on antigen-specific memory T cells in the absence of acute infection could be surprising. However, depending on their specificity, antigen-specific memory T cells can display highly variable differentiation profile [34] and the maintenance of elevated PD-1 expression on some virus-specific memory T cells (such as CMV or HCV) has been reported even after infection clearance [35-37]. Moreover, a recent publication showed that PD-1-expressing SARS-CoV-2-specific CD8 T cells are not exhausted, but functional in convalescent COVID-19 patients [38]. Further experiments will be needed to understand the impact of elevated PD-1 expression of SARS-CoV-2-responding CD4 T cells.

Lastly, SARS-CoV-2-responding CD4 T cells were qualitatively different between convalescent COVID-19 participants and participants who did not experience COVID-19-associated symptoms. In the former group, SARS-CoV-2-specific CD4 T cells displayed almost exclusively an early differentiated memory phenotype; while in the latter group, SARS-CoV-2-responsive CD4 T cells exhibited a late differentiated memory phenotype and were enriched in GrB.

Our study has several limitations. Our results rely on a small number of participants and need to be confirmed in larger cohorts. In the absence of a low SARS-CoV-2 exposure control group, we cannot speculate about the exact specificity of SARS-CoV-2 responding CD4 T cells detected in participants who are serology negative.

Moreover, technical optimization and validation are warranted to streamline the proposed assay and test its performance. A comparable approach using a whole blood IFN- γ release assay (IGRA), commonly used for the detection of *Mycobacterium tuberculosis* infection, has been recently applied to detect SARS-CoV-2 specific T cell response [39]. While both techniques have a comparable cost (~\$30-40), the flow cytometry -based assay proposed in this study has a lower turnaround time compared to IGRA (~7h vs 18h, respectively) but more importantly our approach is more flexible. It can be used to simply quantify SARS-CoV-2 responding CD4 and CD8 T cells (using a minimal 4-color panel) to monitor the frequency of T cell responses in epidemiological studies or vaccine trials or be used for more exploratory work assessing in depth the phenotype and functionality of SARS-CoV-2-specific T cells. Regarding the sensitivity and specificity of each assay, parallel analyses will need to be performed to compare the performance of both assays. Additionally, further optimizations should be performed to improve the standardization and flexibility of the proposed SARS-CoV-2-specific whole blood assay including using tubes containing lyophilized reagents for cell stimulation, pre-mixed antibody cocktail for cell staining and automated gating for flow analyses. Another alternative to alleviate collection time and transport constraints and for limited-resource sites would be to directly cryopreserve whole blood samples using cryopreservative solution, such as CryoStor[®] CS10 solution [40], allowing to perform the cell stimulation and staining in batch at adequately equipped sites. Lastly, in this study, cell stimulation was performed using a pool of peptides covering the SARS-CoV-2 S, N and M proteins. In light of a recent publication showing that the discrimination between pre-existing cross-reactive response and COVID-19-induced response is dependent on specificity of the SARS-CoV-2 antigen used [41], it appears that the peptide combination used in this study is not optimal to distinguish pre-existing cross-reactive SARS-CoV-2-reactive T cells from COVID-19-induced responses. Thus, the proposed assay may be further enhanced by using more targeted peptide pools with defined specificity for SARS-CoV-2 or seasonal coronaviruses.

Nevertheless, the primary objective of this study was to describe, as a proof-of-concept, a rapid and flexible whole blood assay, requiring very limited amount of blood, that could represent a powerful tool to investigate SARS-CoV-2 specific T cell responses both quantitatively and qualitatively even in COVID-19 paediatric cases.

Acknowledgment: The authors thank all the participants.

Funding sources: This work was supported by Wellcome [grant number: 104803 and 203135] and Crick idea to innovation (i2i) scheme in partnership with the Rosetrees Trust [grant number: 2020-0009]. CR is supported by the European and Developing Countries Clinical Trials Partnership EDCTP2 programme supported by the European Union (EU)'s Horizon 2020 programme [grant number: Training and Mobility Action TMA2017SF-1951-TB-SPEC to CR] and the National Institutes of Health (NIH) [grant number:R21AI148027 to CR]. GS is supported by the European and Developing Countries Clinical Trials Partnership EDCTP2 programme [grant number: Training and Mobility Action TMA2018SF-2446 to GS]. RJW and KAW receive support from the Francis Crick institute, which is funded by the UK Medical Research Council UKRI, Cancer Research UK and Wellcome [grant number: FC0010218]. For the purpose of Open Access, the author has applied a CC BY public copyright licence to any Author Accepted Manuscript version arising from this submission.

Author contributions: CR, KAW and RJW designed the study. CR performed the whole blood assay and flow cytometry experiments and analysed and interpreted the data. GS performed the neutralizing assay experiments. EdB, RTG and CS recruited the study participants. KAW, HM and DH provided critical reagents. CR and KAW wrote the manuscript with all authors contributing to providing critical feedback.

Conflict of interest: The authors declare no conflict of interest.

REFERENCES

1. Margolin E, Burgers WA, Sturrock ED, Mendelson M, Chapman R, Douglass N, Williamson AL, Rybicki EP. Prospects for SARS-CoV-2 diagnostics, therapeutics and vaccines in Africa. *Nat Rev Microbiol* 2020: Sep 10;1-15; DOI: 10.1038/s41579-020-00441-3.
2. Cox RJ, Brokstad KA. Not just antibodies: B cells and T cells mediate immunity to COVID-19. *Nat Rev Immunol* 2020: 20(10): 581-582.
3. Altmann DM, Boyton RJ. SARS-CoV-2 T cell immunity: Specificity, function, durability, and role in protection. *Sci Immunol* 2020: 5. DOI: 10.1126/sciimmunol.abd6160(49).
4. Quinti I, Lougaris V, Milito C, Cinetto F, Pecoraro A, Mezzaroma I, Mastroianni CM, Turriziani O, Bondioni MP, Filippini M, Soresina A, Spadaro G, Agostini C, Carsetti R, Plebani A. A possible role for B cells in COVID-19? Lesson from patients with agammaglobulinemia. *J Allergy Clin Immunol* 2020: 146(1): 211-213 e214.
5. Canete PF, Vinuesa CG. COVID-19 makes B cells forget, but T cells remember. *Cell* 2020: 183(1): 13-15.
6. Sette A, Crotty S. Pre-existing immunity to SARS-CoV-2: the knowns and unknowns. *Nat Rev Immunol* 2020: 20(8): 457-458.
7. Lipsitch M, Grad YH, Sette A, Crotty S. Cross-reactive memory T cells and herd immunity to SARS-CoV-2. *Nat Rev Immunol* 2020: 20(11): 709-713.
8. Braun J, Loyal L, Frentsch M, Wendisch D, Georg P, Kurth F, Hippenstiel S, Dingeldey M, Kruse B, Fauchere F, Baysal E, Mangold M, Henze L, Lauster R, Mall MA, Beyer K, Rohmel J, Voigt S, Schmitz J, Miltenyi S, Demuth I, Muller MA, Hocke A, Witzentrath M, Suttorp N, Kern F, Reimer U, Wenschuh H, Drosten C, Corman VM, Giesecke-Thiel C, Sander LE, Thiel A. SARS-CoV-2-reactive T cells in healthy donors and patients with COVID-19. *Nature* 2020: 587(7833):270-274. DOI: 10.1038/s41586-020-2598-9.
9. Mendelson M, Boloko L, Boutall A, Cairncross L, Calligaro G, Coccia C, Dave JA, de Villiers M, Dlamini S, Frankenfeld P, Gina P, Gule MV, Hoare J, Hofmeyr R, Hsiao M, Joubert I, Kahn T, Krause R, Kroopman A, Levin D, Maughan D, Mazondwa S, Meintjes G, Nordien R, Ntusi N, Papavarnavas N, Peter J, Pickard H, Raubenheimer P, Said-Hartley Q, Singh P, Wasserman S, Team obotGSHC-R. Clinical management of COVID-19: Experiences of the COVID-19 epidemic from Groote Schuur Hospital, Cape Town, South Africa. *South African Medical Journal* 2020: 110: 973-981.
10. Muench P, Jochum S, Wenderoth V, Ofenloch-Haehnle B, Hombach M, Strobl M, Sadlowski H, Sachse C, Torriani G, Eckerle I, Riedel A. Development and validation of the Elecsys anti-SARS-CoV-2 immunoassay as a highly specific tool for determining past exposure to SARS-CoV-2. *J Clin Microbiol* 2020: 58. DOI: 10.1128/JCM.01694-20(10).
11. Favresse J, Eucher C, Elsen M, Tre-Hardy M, Dogne JM, Douxfils J. Clinical performance of the Elecsys electrochemiluminescent immunoassay for the detection of SARS-CoV-2 total antibodies. *Clin Chem* 2020: 66(8): 1104-1106.
12. National S-C-SAEG. Performance characteristics of five immunoassays for SARS-CoV-2: a head-to-head benchmark comparison. *Lancet Infect Dis* 2020: 20(12):1390-1400 DOI: 10.1016/S1473-3099(20)30634-4.
13. Buck CB, Pastrana DV, Lowy DR, Schiller JT. Efficient intracellular assembly of papillomaviral vectors. *J Virol* 2004: 78(2): 751-757.
14. Buck CB, Pastrana DV, Lowy DR, Schiller JT. Generation of HPV pseudovirions using transfection and their use in neutralization assays. *Methods Mol Med* 2005: 119: 445-462.
15. Rogers TF, Zhao F, Huang D, Beutler N, Burns A, He WT, Limbo O, Smith C, Song G, Woehl J, Yang L, Abbott RK, Callaghan S, Garcia E, Hurtado J, Parren M, Peng L, Ramirez S, Ricketts J, Ricciardi MJ, Rawlings SA, Wu NC, Yuan M, Smith DM, Nemazee D, Tejjaro JR, Voss JE, Wilson IA, Andrabi R, Briney B, Landais E, Sok D, Jardine JG, Burton DR.

- Isolation of potent SARS-CoV-2 neutralizing antibodies and protection from disease in a small animal model. *Science* 2020: 369(6506): 956-963.
16. Mou H, Quinlan BD, Peng H, Guo Y, Peng S, Zhang L, Davis-Gardner ME, Gardner MR, Crynen G, Voo ZX, Bailey CC, Alpert MD, Rader C, Choe H, Farzan M. Mutations from bat ACE2 orthologs markedly enhance ACE2-Fc neutralization of SARS-CoV-2. *BioRxiv* 2020(10.1101/2020.06.29.178459): DOI: 10.1101/2020.1106.1129.178459.
 17. Hanekom WA, Hughes J, Mavinkurve M, Mendillo M, Watkins M, Gamieldien H, Gelderbloem SJ, Sidibana M, Mansoor N, Davids V, Murray RA, Hawkridge A, Haslett PA, Ress S, Hussey GD, Kaplan G. Novel application of a whole blood intracellular cytokine detection assay to quantitate specific T-cell frequency in field studies. *J Immunol Methods* 2004: 291(1-2): 185-195.
 18. Riou C, Du Bruyn E, Ruzive S, Goliath RT, Lindestam Arlehamn CS, Sette A, Sher A, Barber DL, Wilkinson RJ. Disease extent and anti-tubercular treatment response correlates with *Mycobacterium tuberculosis*-specific CD4 T-cell phenotype regardless of HIV-1 status. *Clin Transl Immunology* 2020: 9(9): e1176.
 19. Liu Y, Yan LM, Wan L, Xiang TX, Le A, Liu JM, Peiris M, Poon LLM, Zhang W. Viral dynamics in mild and severe cases of COVID-19. *Lancet Infect Dis* 2020: 20(6): 656-657.
 20. Liu L, Liu W, Zheng Y, Jiang X, Kou G, Ding J, Wang Q, Huang Q, Ding Y, Ni W, Wu W, Tang S, Tan L, Hu Z, Xu W, Zhang Y, Zhang B, Tang Z, Zhang X, Li H, Rao Z, Jiang H, Ren X, Wang S, Zheng S. A preliminary study on serological assay for severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) in 238 admitted hospital patients. *Microbes Infect* 2020: 22(4-5): 206-211.
 21. Wang J, Chen C, Li Q, Cai P, Wang Z, Wang L. COVID-19 confirmed patients with negative antibodies results. *BMC Infect Dis* 2020: 20(1): 698.
 22. Juno JA, van Bockel D, Kent SJ, Kelleher AD, Zaunders JJ, Munier CM. Cytotoxic CD4 T cells - Friend or foe during viral infection? *Front Immunol* 2017: 8: 19.
 23. Pachnio A, Ciauriz M, Begum J, Lal N, Zuo J, Beggs A, Moss P. Cytomegalovirus infection leads to development of high frequencies of cytotoxic virus-specific CD4+ T cells targeted to vascular endothelium. *PLoS Pathog* 2016: 12(9): e1005832.
 24. Mathew D, Giles JR, Baxter AE, Oldridge DA, Greenplate AR, Wu JE, Alanio C, Kuri-Cervantes L, Pampena MB, D'Andrea K, Manne S, Chen Z, Huang YJ, Reilly JP, Weisman AR, Ittner CAG, Kuthuru O, Dougherty J, Nzingha K, Han N, Kim J, Pattekar A, Goodwin EC, Anderson EM, Weirick ME, Gouma S, Arevalo CP, Bolton MJ, Chen F, Lacey SF, Ramage H, Cherry S, Hensley SE, Apostolidis SA, Huang AC, Vella LA, Unit UCP, Betts MR, Meyer NJ, Wherry EJ. Deep immune profiling of COVID-19 patients reveals distinct immunotypes with therapeutic implications. *Science* 2020: 369(6508):eabc8511. DOI: 10.1126/science.abc8511.
 25. Chen Z, Wherry EJ. T cell responses in patients with COVID-19. *Nat Rev Immunol* 2020: 20(9): 529-536.
 26. Grifoni A, Weiskopf D, Ramirez SI, Mateus J, Dan JM, Moderbacher CR, Rawlings SA, Sutherland A, Premkumar L, Jadi RS, Marrama D, de Silva AM, Frazier A, Carlin AF, Greenbaum JA, Peters B, Krammer F, Smith DM, Crotty S, Sette A. Targets of T cell responses to SARS-CoV-2 coronavirus in humans with COVID-19 disease and unexposed individuals. *Cell* 2020: 181(7): 1489-1501 e1415.
 27. Sekine T, Perez-Potti A, Rivera-Ballesteros O, Stralin K, Gorin JB, Olsson A, Llewellyn-Lacey S, Kamal H, Bogdanovic G, Muschiol S, Wullimann DJ, Kammann T, Emgard J, Parrot T, Folkesson E, Karolinska C-SG, Rooyackers O, Eriksson LI, Henter JI, Sonnerborg A, Allander T, Albert J, Nielsen M, Klingstrom J, Gredmark-Russ S, Bjorkstrom NK, Sandberg JK, Price DA, Ljunggren HG, Aleman S, Buggert M. Robust T cell immunity in convalescent individuals with asymptomatic or mild COVID-19. *Cell* 2020: 183(1): 158-168 e114.

28. Le Bert N, Tan AT, Kunasegaran K, Tham CYL, Hafezi M, Chia A, Chng MHY, Lin M, Tan N, Linster M, Chia WN, Chen MI, Wang LF, Ooi EE, Kalimuddin S, Tambyah PA, Low JG, Tan YJ, Bertoletti A. SARS-CoV-2-specific T cell immunity in cases of COVID-19 and SARS, and uninfected controls. *Nature* 2020: 584(7821): 457-462.
29. Rydyznski Moderbacher C, Ramirez SI, Dan JM, Grifoni A, Hastie KM, Weiskopf D, Belanger S, Abbott RK, Kim C, Choi J, Kato Y, Crotty EG, Kim C, Rawlings SA, Mateus J, Tse LPV, Frazier A, Baric R, Peters B, Greenbaum J, Ollmann Saphire E, Smith DM, Sette A, Crotty S. Antigen-specific adaptive immunity to SARS-CoV-2 in acute COVID-19 and associations with age and disease severity. *Cell* 2020: 183(4):996-1012.e19. DOI: 10.1016/j.cell.2020.09.038.
30. Blanco-Melo D, Nilsson-Payant BE, Liu WC, Uhl S, Hoagland D, Moller R, Jordan TX, Oishi K, Panis M, Sachs D, Wang TT, Schwartz RE, Lim JK, Albrecht RA, tenOever BR. Imbalanced host response to SARS-CoV-2 drives development of COVID-19. *Cell* 2020: 181(5): 1036-1045 e1039.
31. Sattler A, Angermair S, Stockmann H, Heim KM, Khadzhynov D, Treskatsch S, Halleck F, Kreis ME, Kotsch K. SARS-CoV-2 specific T-cell responses and correlations with COVID-19 patient predisposition. *J Clin Invest* 2020: 130(12):6477-6489. DOI: 10.1172/JCI140965.
32. Schub D, Klemis V, Schneitler S, Mihm J, Lepper PM, Wilkens H, Bals R, Eichler H, Gartner BC, Becker SL, Sester U, Sester M, Schmidt T. High levels of SARS-CoV-2-specific T cells with restricted functionality in severe courses of COVID-19. *JCI Insight* 2020: 5(20): e142167.
33. Zheng HY, Zhang M, Yang CX, Zhang N, Wang XC, Yang XP, Dong XQ, Zheng YT. Elevated exhaustion levels and reduced functional diversity of T cells in peripheral blood may predict severe progression in COVID-19 patients. *Cell Mol Immunol* 2020: 17(5): 541-543.
34. Appay V, Dunbar PR, Callan M, Klenerman P, Gillespie GM, Papagno L, Ogg GS, King A, Lechner F, Spina CA, Little S, Havlir DV, Richman DD, Gruener N, Pape G, Waters A, Easterbrook P, Salio M, Cerundolo V, McMichael AJ, Rowland-Jones SL. Memory CD8+ T cells vary in differentiation phenotype in different persistent virus infections. *Nat Med* 2002: 8(4): 379-385.
35. Parry HM, Dowell AC, Zuo J, Verma K, Kinsella FAM, Begum J, Croft W, Sharma-Oates A, Pratt G, Moss P. PD-1 is imprinted on cytomegalovirus-specific CD4+ T cells and attenuates Th1 cytokine production whilst maintaining cytotoxicity. *PLoS Pathog* 2021: 17(3): e1009349.
36. Kasprowicz V, Schulze Zur Wiesch J, Kuntzen T, Nolan BE, Longworth S, Berical A, Blum J, McMahan C, Reyor LL, Elias N, Kwok WW, McGovern BG, Freeman G, Chung RT, Klenerman P, Lewis-Ximenez L, Walker BD, Allen TM, Kim AY, Lauer GM. High level of PD-1 expression on hepatitis C virus (HCV)-specific CD8+ and CD4+ T cells during acute HCV infection, irrespective of clinical outcome. *J Virol* 2008: 82(6): 3154-3160.
37. Jubel JM, Barbati ZR, Burger C, Wirtz DC, Schildberg FA. The role of PD-1 in acute and chronic infection. *Front Immunol* 2020: 11: 487.
38. Rha MS, Jeong HW, Ko JH, Choi SJ, Seo IH, Lee JS, Sa M, Kim AR, Joo EJ, Ahn JY, Kim JH, Song KH, Kim ES, Oh DH, Ahn MY, Choi HK, Jeon JH, Choi JP, Kim HB, Kim YK, Park SH, Choi WS, Choi JY, Peck KR, Shin EC. PD-1-Expressing SARS-CoV-2-Specific CD8(+) T Cells Are Not Exhausted, but Functional in Patients with COVID-19. *Immunity* 2021: 54(1): 44-52 e43.
39. Murugesan K, Jagannathan P, Pham TD, Pandey S, Bonilla HF, Jacobson K, Parsonnet J, Andrews JR, Weiskopf D, Sette A, Pinsky BA, Singh U, Banaei N. Interferon-gamma release assay for accurate detection of SARS-CoV-2 T cell response. *Clin Infect Dis* 2020. Oct 9;ciaa1537. doi: 10.1093/cid/ciaa1537.

40. Braudeau C, Salabert-Le Guen N, Chevreuil J, Rimbart M, Martin JC, Josien R. An easy and reliable whole blood freezing method for flow cytometry immuno-phenotyping and functional analyses. *Cytometry B Clin Cytom* 2021, Feb 5. doi: 10.1002/cyto.b.21994.
41. Ogbe A, Kronsteiner B, Skelly DT, Pace M, Brown A, Adland E, Adair K, Akhter HD, Ali M, Ali SE, Angyal A, Ansari MA, Arancibia-Carcamo CV, Brown H, Chinnakannan S, Conlon C, de Lara C, de Silva T, Dold C, Dong T, Donnison T, Eyre D, Flaxman A, Fletcher H, Gardner J, Grist JT, Hackstein CP, Jaruthamsophon K, Jeffery K, Lambe T, Lee L, Li W, Lim N, Matthews PC, Mentzer AJ, Moore SC, Naisbitt DJ, Ogese M, Ogg G, Openshaw P, Pirmohamed M, Pollard AJ, Ramamurthy N, Rongkard P, Rowland-Jones S, Sampson O, Screatton G, Sette A, Stafford L, Thompson C, Thomson PJ, Thwaites R, Vieira V, Weiskopf D, Zacharopoulou P, Oxford Immunology Network Covid-19 Response TCC, Oxford Protective TCIc-CT, Turtle L, Klenerman P, Goulder P, Frater J, Barnes E, Dunachie S. T cell assays differentiate clinical and subclinical SARS-CoV-2 infections from cross-reactive antiviral responses. *Nat Commun* 2021; 12(1): 2055.

FIGURE LEGENDS

Figure 1: Schematic showing methodology and workflow of the whole blood assay for the detection of SARS-CoV-2-specific adaptive immune responses. Step 1: 400 µl of heparinized whole blood is incubated for 5 hours in the presence of a SARS-CoV-2-specific peptide pool in the presence of co-stimulatory antibodies (i.e. CD28 and CD49d) and Brefeldin-A. **Step 2:** Cells are incubated for 20 min in the presence of a transcription factor fixation buffer, leading to the simultaneous lysis of red blood cells and cell fixation. **Step 3:** Cells are stained for 30 min with an optimized panel of fluorophore labelled antibodies. **Step 4:** Samples are acquired on a flow cytometer. Control samples are processed with a similar workflow in the absence of SARS-CoV-2-specific peptide pool.

Figure 1




Figure 2: SARS-CoV-2 serological assessment. A- Quantification of SARS-CoV-2 Nucleocapsid-specific antibodies using the Elecsys[®] Roche assay expressed as a cut-off index (signal sample/cut-off). Participants were grouped according to their clinical characteristics (Blue: no symptoms, no SARS-CoV-2 PCR performed; Orange: self-reported symptoms, SARS-CoV-2 PCR negative; and Red: self-reported symptoms, SARS-CoV-2 PCR positive). Bars represent the medians. The dotted line indicates the manufacturer's cut-off value for positivity. Statistical comparisons were performed using a Mann-Whitney T- test. **B-** SARS-CoV-2 pseudovirus neutralization activity. SARS-CoV-2 pseudovirions pre-incubated with serially diluted patient plasma were used to infect ACE2-expressing HEK-293T cells. Luciferase activity as a measure for infection was assessed 3 days post-infection, and results are expressed as infection compared to control (untreated virions, grey shaded area) which was set at 100%.

Figure 2




Figure 3: Magnitude and functional profile of SARS-CoV-2-specific CD4 T cells. A-

Proportion of participants exhibiting a detectable SARS-CoV-2-specific CD4 T cell response. **B-** Representative examples of TNF- α and IFN- γ production in CD4 T cells in response to SARS-CoV-2 peptide pool. NS: no stimulation. **C-** Magnitude of SARS-CoV-2-specific CD4 T cell response (expressed as a percentage of total CD4 T cells) in participants grouped according to their clinical characteristics (Blue: no symptoms, no SARS-CoV-2 PCR performed; Orange: self-reported symptoms, SARS-CoV-2 PCR negative; and Red: self-reported symptoms, SARS-CoV-2 PCR positive). The number of participants and % of responders in each group is presented at the bottom of the graph. Statistical comparisons were performed using the Kruskal-Wallis test. **D-** Polyfunctional profile of SARS-CoV-2-specific CD4 T cells in each group. The x-axis displays the composition of each combination which is denoted with a dot for the presence of IL-2, IFN- γ and TNF- α . The medians (black bar) are shown. Each combination is color-coded, and data are summarized in the pie charts, where each pie slice represents the median contribution of each combination to the total SARS-CoV-2 response. The arcs identify the contribution of TNF α (light grey), IL-2 (grey) and IFN- γ (black) to SARS-CoV-2 response. The Wilcoxon Rank Sum Test was used to compare response patterns between groups. Statistical differences between pie charts were defined using a permutation test.




Figure 4: Memory differentiation profile and Granzyme B (GrB) expression in SARS-CoV-2-specific CD4 T cells. **A-** Representative examples of the memory differentiation profile of SARS-CoV-2-specific CD4 T cells based on the expression of CD45RA and CD27. The flow plot on the right shows the distribution of Naïve (CD45RA+CD27+), Early differentiated (ED: CD45RA-CD27+), Late differentiated (LD: CD45RA-CD27-) and Effector (Eff: CD45RA+CD27-) in total CD4 T cells. **B-** Summary graph of the proportion of ED and LD in SARS-CoV-2-specific CD4 T cells in each group (Blue: no symptoms, no SARS-CoV-2 PCR performed; Orange: self-reported symptoms, SARS-CoV-2 PCR negative and Red: self-reported symptoms, SARS-CoV-2 PCR positive). Statistical comparisons were performed using the Kruskal-Wallis test. **C-** Representative examples of GrB expression in total and SARS-CoV-2-specific CD4 T cells. **D-** Summary graph of GrB expression in SARS-CoV-2-specific CD4 T cells in each group. Statistical comparisons were performed using the Kruskal-Wallis test. **E-** Relationship between the proportion of ED within SARS-CoV-2-specific CD4 T cells and GrB expression or the proportion of IFN γ +TNF α +IL2- SARS-CoV-2-specific CD4 T cells. Correlations were tested by a two-tailed non-parametric Spearman rank test.

Figure 4




Figure 5: Phenotypic signature of SARS-CoV-2-specific IFN γ + CD4+ T cells according to clinical characteristics. **A-** Non-supervised two-way hierarchical cluster analysis (HCA, Ward method) using three phenotypic parameters (i.e. GrB expression and the proportion of ED and IFN γ +TNF α +IL2-) from SARS-CoV-2-specific CD4 T cells. Each column represents a participant and is color-coded according to their clinical characteristics indicated by a dot at the top of the dendrogram (Blue: no symptoms, no SARS-CoV-2 PCR performed; Orange: self-reported symptoms, SARS-CoV-2 PCR negative; and Red: self-reported symptoms, SARS-CoV-2 PCR positive). Participants with a positive SARS-CoV-2 serology test are indicated by a black box and a white box identifies SARS-CoV-2 antibody-negative subjects. Data are depicted as a heatmap coloured from minimum to maximum values for each parameter. **B-** Principal component analysis on correlations, derived from the three studied parameters. Each dot represents a participant. The two axes represent principal components 1 (PC1) and 2 (PC2). Their contribution to the total data variance is shown as a percentage. **C-** Loading plot showing how each parameter influences PC1 and PC2 values.







Figure 6: Activation profile of SARS-CoV-2-specific CD4 T cells. A- CD38, HLA-DR, Ki67 and PD-1 expression in SARS-CoV-2-specific CD4 T cells in each group (Blue: no symptoms, no SARS-CoV-2 PCR performed; Orange: self-reported symptoms, SARS-CoV-2 PCR negative and Red: self-reported symptoms, SARS-CoV-2 PCR positive). No significant differences were observed between groups for any markers, using the Kruskal-Wallis test. **B-** Association between CD38 expression in SARS-CoV-2-specific CD4 T cells and the time post SARS-CoV-2 PCR positive test (weeks). Each symbol represents a participant (n=9). Dotted red lines identify participant with longitudinal samples. Correlations were tested by a two-tailed non-parametric Spearman rank test. **C-** Comparison of the correlation between the time post SARS-CoV-2 PCR positive test (weeks) and the expression of different activation profile markers, ranked according to the strength of the association. Spearman correlation r values are plotted on the x-axis and corresponding P -values are shown within each bar.



	Healthcare workers (n=31)			<i>P</i> -values
n	15	7	9	na
Self-reported COVID-19-like symptoms	No	Yes	Yes	na
Tested for COVID-19 (PCR)	No	Yes	Yes	na
Positive SARS-CoV-2 PCR	na	No	Yes	na
Positive SARS-CoV-2 N-spe IgG (%) ^a	0%	0%	88.9% (8/9)	na
Presence of SARS-CoV-2 nAbs (%) ^b	0%	0%	100%	na
Contact with COVID-19 patients (%)	86.6%	100%	88.9%	NS
Time post PCR test (weeks) ^c	na	7.3 [6.3-8.4]	4.7 [2.6-7.2]	0.09

Table 1. Clinical characteristics of study participants

^a Nucleocapsid-specific IgG measured using the Elecsys Roche system.

^b SARS-CoV-2 neutralizing antibodies (nAbs) measured using a SARS-CoV-2 pseudovirus neutralisation assay.

^c median and interquartile range (IQR)

na: not applicable, NS: non-significant

Online supplement for:

**Rapid, simplified whole blood-based multiparameter assay to quantify and phenotype
SARS-CoV-2 specific T cells**

Catherine Riou^{1,2,*}, Georgia Schäfer^{1,3,4}, Elsa du Bruyn^{1,5}, Rene T. Goliath¹, Cari Stek^{1,5,6}, Huihui Mou⁷, Deli Hung⁸, Katalin A. Wilkinson^{1,5,9}, Robert J. Wilkinson^{1,5,6,9}

¹ Wellcome Centre for Infectious Disease Research in Africa and Institute of Infectious Disease and Molecular Medicine, University of Cape Town, Observatory 7925, South Africa

² Division of Immunology, Department of Pathology, University of Cape Town, Observatory 7925, South Africa

³ International Centre for Genetic Engineering and Biotechnology (ICGEB) Cape Town, Observatory 7925, South Africa

⁴ Division of Medical Biochemistry and Structural Biology, Department of Integrative Biomedical Sciences, University of Cape Town, Observatory 7925, South Africa

⁵ Department of Medicine, University of Cape Town, Observatory 7925, South Africa

⁶ Department of Infectious Diseases, Imperial College London, W12 0NN, UK


⁷ Department of Immunology and Microbiology, The Scripps Research Institute, Jupiter, FL 33458, USA

⁸ Department of Immunology and Microbiology, The Scripps Research Institute, La Jolla, CA 92037, USA


⁹ The Francis Crick Institute, London, NW1 1AT, UK

Markers	Fluorochrome	Clone	Company	Cat. Number	Role
CD3	BV650	OKT3	BioLegend	317323	Lineage
CD4	BV785	OKT4	BioLegend	317428	
CD8	BV510	RPA-T8	BioLegend	301048	
CD45RA	Alexa 488	HI100	BioLegend	304114	Memory differentiation
CD27	PE-Cy5	1A4CD27	Beckman	6607107	
CD38	APC	HIT2	BD Bioscience	555462	Activation
HLA-DR	BV605	L243	BioLegend	307640	
Ki67	PerCP-Cy5.5	B56	BD Bioscience	561284	
PD-1	PE	EH12.2H7	BioLegend	329906	
GrB	BV421	GB11	BD Bioscience	563388	Cytotoxic potential
IFNγ	BV711	4S.B3	BioLegend	502540	Functions
TNFα	PE-Cy7	MAB11	BioLegend	502930	
IL-2	PE/Dazzle™ 594	MQ1-17H12	BioLegend	500344	

Supp table 1: Description of the antibody panel used in the study.



Supp Figure 1: Gating strategy used to identify SARS-CoV-2 specific CD4 T cells and defined their phenotypic characteristics. **A-** Cell lineage. **B-** Phenotyping of total CD4 T cell and cytokine production in response to SARS-CoV-2 peptides. **C-** Phenotype of SARS-CoV-2 responding CD4 T cells in one participant with no symptom (left) and one COVID-19 confirmed participant (right). SARS-CoV-2-specific CD4 T cells (expressing any measured cytokine) are depicted in red.



Supp Figure 2: Comparison of magnitude of SARS-CoV-2-specific CD4 T cell response expressing TNF α , IL-2 or IFN γ in participants exhibiting a detectable response to SARS-CoV-2 peptide pool. Each dot represent a participant and is color-coded according to clinical characteristics (Blue represents no symptoms, no SARS-CoV-2 PCR performed; Orange: self-reported symptoms, SARS-CoV-2 PCR negative and Red: self-reported symptoms, SARS-CoV-2 PCR positive). Black bars indicate medians. Statistical comparisons were performed using a nonparametric paired Friedman test.